










(IT) **Istruzioni di lavaggio:** Detergenti neutri e non abrasivi - Autoclave a 85 °C per un minuto - (EN) **Washing instructions:** Detergent neutral and non-abrasive - autoclave at 85 °C for one minute -
(FR) **Instructions de lavage:** Détergent neutre et non abrasif - autoclave à 85 °C pendant une minute - (E) **Instrucciones de lavado:** Detergente neutro y no abrasivo - autoclave a 85 °C durante un minuto - (D) **Waschanleitung:** Waschmittel neutral und nicht-abrasive - Autoklaven bei 85 °C für eine Minute - (PT) **Instruções de lavagem:** Detergente neutro e não abrasivo - autoclavável a 85°C durante 1 minuto.